

CANI, Inc. Detergent Selection Chart

Detergent Type	Cleaning Methods			Cleaning Agent Cross Reference
	Manual, Immersion	Machine Washers	Recirculating CIP	
Alkaline	X	X	X	CIP 100®
CP-700				
Chlorinated Alkaline				
CP-725	X	X	X	CIP 150
Acidic Inorganic				
CP-450	X	X	X	CIP 200
Acidic Organic				
CP-400	X	X	X	CIP 220
Neutral pH				
CP-830	X	X	X	CIP 300
Microbial Control				
Minncare	X	X	X	Spor-Klenz
Actril	X			Spor-Klenz RTU